

0.5% SERUM Kdo2 LIPID A TREATMENT PROTOCOL
Version 1

1. See Figure 1 for an overview of the time course of the protocol.
2. Maintain sterile technique throughout the protocol until harvesting.
3. Plate cells as recommended in the LIPID MAPS Thawing and Passage Procedure in growth medium to achieve 80% confluence at 30 hours following plating.
4. Assign a barcode to each plate/vessel and enter into LIMS.
5. Make up Serum Deprivation Media (DMEM, 0.5% LIPID MAPS FBS; LIPID MAPS solution protocol PS0000001000). Assign a barcode to each bottle and enter into LIMS.
6. Thirty hours after plating, rinse the vessels 1 X with 37° C DPBS.
7. Add a volume of fresh 37° C Serum Deprivation Medium.
8. Incubate for 18 hours at 37° C.
9. Spray the Eppendorf containing the freshly sonicated 1000x working solution (100 µg/ml) of Kdo2 Lipid A with 70% ethanol and let air dry before using.
10. Add the Kdo2 Lipid A to the medium, for a final concentration of 100 ng/ml. Add an equal amount of DPBS to controls.
11. Incubate for 24 hours at 37° C.
12. Before harvesting the cells, take an aliquot (1 ml) of medium from each condition for the TNF α assay. Assign a barcode and enter into LIMS. Follow the TNF α assay protocol in R&D's manual.
For the TNF α assay;
The TNF α aliquots from the LPS induced cells must be diluted in medium 1:20, 1:40 and 1:80 before assaying. Do not dilute aliquots from cells that were not treated with Kdo2 lipid A.
13. Place vessels on ice, aspirate, and wash each vessel 2 X with an appropriate volume of 4° C DPBS.
14. Add another fresh volume of 4° C DPBS and scrape the cells with a scraper (see equipment list).
15. Pipette the cell suspension into an appropriate tube for either direct lipid extraction or centrifugation. Assign a barcode and enter into LIMS.
16. Suspend the cells and take an aliquot, e.g., 200 µl for 20 µl duplicates, for DNA analysis. Assign a barcode and enter into LIMS. Aliquots can be frozen for later DNA analysis. Follow the DNA assay protocol in Molecular Probe's manual..
17. Cells can now be extracted directly or spun down for extraction of cell pellets. To centrifuge cells, spin the cell suspension at 2000 rpm for 5-10 minutes at 4° C.

Figure 1

