

**10% SERUM Kdo<sub>2</sub> LIPID A TREATMENT PROTOCOL**  
**Version 3, 8-26-05**

1. See Figure 1 for an overview of the time course of the Kdo<sub>2</sub> induction protocol.
2. Maintain sterile technique throughout the Kdo<sub>2</sub> induction protocol until harvesting.
3. Plate cells as recommended on the LIPID MAPS Thawing and Passage Procedure in growth medium to achieve 80% confluence at 30 hours following plating.
4. Assign a barcode to each plate/vessel and enter into LIMS.
5. Thirty hours after plating, rinse the vessels 1 X with 37° C DPBS.
6. Add a volume of fresh 37° C growth medium.
7. Incubate for 18 hours at 37° C.
8. Immediately after removing the plate from 37° C and before treating the cells, take an aliquot (1 ml) of medium from each condition for the TNF $\alpha$  assay. Assign a barcode, enter into LIMS and freeze the aliquots at -20° C.
9. Spray the Eppendorf containing the freshly sonicated 1000x working solution (100  $\mu$ g/ml) of Kdo<sub>2</sub> lipid A with 70% ethanol and let air dry before using.
10. Add the Kdo<sub>2</sub> Lipid A to the medium, for a final concentration of 100 ng/ml. Add an equal amount of DPBS to controls.
11. Incubate for 24 hours at 37° C.
12. Immediately after removing the plate from 37° C and before harvesting the cells, take an aliquot (1 ml) of medium from each condition for the TNF $\alpha$  assay. Assign a barcode, enter into LIMS and freeze the aliquots at -20° C.  
For the TNF $\alpha$  assay;  
The TNF $\alpha$  aliquots from the Kdo<sub>2</sub> induced cells must be diluted in medium 1:40 and 1:80 before assaying. Do not dilute aliquots from cells that were not treated with Kdo<sub>2</sub> lipid A. Send frozen aliquots to ElisaTech, 12635 E. Montview Blvd., Suite 216, Aurora, CO 80010.
13. Place vessels on ice, aspirate, and wash each vessel 2 X with an appropriate volume of 4° C DPBS.
14. Add another fresh volume of 4° C DPBS and scrape the cells with a scraper (see equipment list).
15. Pipet the cell suspension into an appropriate tube for either direct lipid extraction or centrifugation. Assign a barcode and enter into LIMS.
16. Suspend the cells and take an aliquot, e.g., 200  $\mu$ l for 20  $\mu$ l duplicates, for DNA analysis. Assign a barcode and enter into LIMS. Aliquots can be frozen for later DNA analysis. Follow the DNA assay protocol in Molecular Probe's manual with the exception of using 5  $\mu$ l of standards instead of 10  $\mu$ l. If you anticipate having a lower DNA concentration, use less standard for your DNA curve.
17. Cells can now be extracted directly or spun down for extraction of cell pellets. To centrifuge cells, spin the cell suspension at 2000 rpm for 5-10 minutes at 4° C.

Figure 1



