

LIPID MAPS Solution Protocol

Solution name: 25-Hydroxycholesterol

Solution name abbreviation: 25OHC

Protocol ID:

Version: 00

Volume:

Components:

Reagent	Source	Catalog or Protocol No.	F.W. or Stock Conc.	Quantity	Final Concentration
25OHC	Sigma	H1015	5mM		5uM
DMSO	Sigma	D2650	100%		

Solution Preparation:

1. Dissolve 25OHC in DMSO.
2. Aliquot in 1.5 ml labeled Eppendorf tubes.
3. Replace air above solution with Argon before freezing at -20°C.

Author: Donna Reichart

Date: 3-22-10